

Aim:

To determine the concentration of the given hydrochloric acid (HCl) solution by titrating it against a standard sodium carbonate (Na_2CO_3) solution.

Apparatus and Chemicals Required:**Apparatus:**

- Burette
- Pipette (10 mL)
- Conical flask
- Beaker
- Funnel
- Stand with clamp
- White tile (optional)

Chemicals:

- Sodium carbonate (Na_2CO_3) solution
 - Hydrochloric acid (HCl) solution
 - Methyl orange indicator
 - Distilled water
-

Procedure:**1. Preparation of Sodium Carbonate Solution:**

- Weigh a known mass of anhydrous sodium carbonate (Na_2CO_3) and dissolve it in distilled water.
- Transfer the solution to a volumetric flask and make up the volume to a known mark.

2. Filling the Burette:

- Rinse the burette with distilled water and then with hydrochloric acid.
- Fill the burette with hydrochloric acid and record the initial reading.

3. Pipetting the Sodium Carbonate Solution:

- Rinse the pipette with distilled water and then with sodium carbonate solution.
- Use a pipette to transfer 10 mL of sodium carbonate solution into a conical flask.
- Add 2-3 drops of methyl orange indicator (solution turns yellow).

4. Performing the Titration:

- Slowly add hydrochloric acid from the burette to the conical flask while swirling it continuously.
- Near the endpoint, add acid dropwise until the color changes from yellow to orange-pink.
- Record the final burette reading.

5. Repeating the Experiment:

- Perform at least three titrations and take the average volume of hydrochloric acid used.

Result:

The concentration of the given hydrochloric acid solution is determined using the titration data.

Balanced Chemical Equation:

